




### CT740FCN-774VR-QR Pressure Regulator Fixture, Fume Hood Mounted

**THIS FIXTURE IS FOR USE ON PURE GAS SYSTEMS, INCLUDING OXYGEN.** THE FOLLOWING GASES MAY BE USED WITH THIS FIXTURE: AIR, ARGON, CARBON DIOXIDE, HELIUM, NITROGEN, NITROUS OXIDE AND OXYGEN. DO NOT USE WHERE PRESSURE OR TEMPERATURE CAN EXCEED RATED OPERATING CONDITIONS.



#### NOTES:

1. VALVE IS FACTORY ASSEMBLED AND TESTED AT 375 PSI HELIUM PRESSURE. MAXIMUM WORKING PRESSURE IS 250 PSI.
2. FIXTURE IS CLEANED, DEGREASED AND PACKED FOR PURE GAS SERVICE.
3. PROVIDE  $\phi 1-1/4"$  MOUNTING HOLE FOR VALVE.
4. PANEL THICKNESS NOT TO EXCEED 1/8".

Measurements may vary  $\pm 1/4"$ .

© 2008 / WaterSaver Faucet Co.