




- **L3100-131WSA** Laboratory Water Valve Assembly, Overhead Mounted Single
- **L3100-110-131WSA** Same as Above Except with L110 Vacuum Breaker




**L3100-131WSA**  
Single




**L3100-132AWSA**  
90° Double





**L3100-132SWSA**  
180° Double



**L3100-133WSA**  
3-Way



**L3100-134WSA**  
4-Way



WATER SAVER RENEWABLE UNIT WITH REPLACEABLE STAINLESS STEEL SEAT



NOTES:

1. FIXTURE IS CERTIFIED BY CSA INTERNATIONAL TO COMPLY WITH ASME A112.18.1M AND CSA B125.1.
2. FIXTURE IS FULLY ASSEMBLED AND FACTORY TESTED PRIOR TO SHIPMENT.

Measurements may vary ± 1/4".

© 2008 / WaterSaver Faucet Co.