

○ L3173-366-158WS Pressure Regulator Fixture, Panel Mounted

THIS FIXTURE IS FOR USE ON PURE GAS SYSTEMS, INCLUDING OXYGEN. THE FOLLOWING GASES MAY BE USED WITH THIS FIXTURE: AIR, ARGON, CARBON DIOXIDE, HELIUM, NITROGEN, NITROUS OXIDE AND OXYGEN. DO NOT USE WHERE PRESSURE OR TEMPERATURE CAN EXCEED RATED OPERATING CONDITIONS.


PRESSURE REGULATOR SPECIFICATIONS:

INLET PRESSURE: 250 PSIG Max.

OUTLET PRESSURE RANGE: 2-110 PSIG

TEMPERATURE RANGE: 0°-125°F

PISTON: NON-RELIEVING PLASTIC
WITH BUNA-N SEALS



Measurements may vary $\pm 1/4$ ".

© 2008 / WaterSaver Faucet Co.